

FR2+ Defibrillator



M3860A WITH ECG and manual override



M3861A WITHOUT ECG or manual override

The Philips HeartStart FR2+ has come to be a trusted tool of professional responders and designated response teams in the work place and public settings. Its clean, uncluttered design is optimized for fast, efficient operation and rapid delivery of defibrillation therapy. Commands are clear and concise: ideal for responders who are trained, drilled, and ready to help save a life now.

The HeartStart FR2+ offers features that make hand-off to advanced life support (ALS) trained professionals seamless to ensure continuity of care for the patient. An ECG display is available**, and the FR2+ can be configured to allow ALS responders to switch to a manual mode**, giving them more decision-making control. And defibrillator pads adapters enable the pads to remain on the patient when transferring to many popular ALS manual defibrillators from Philips and other manufacturers.

But don't be misled by its advanced capabilities. The HeartStart FR2+ is extremely easy to use. Its design is based on years of research refined through Philips human factors experts, and user testing with the full range of responders, from healthcare professionals to lay people. In fact, hundreds of thousands of citizens are already trained and ready to help save a life with the FR2+.

Optimized for challenging environments

The HeartStart FR2+ is the right choice for harsh, outdoor, or mobile use. It combines natural sounding voice instructions that are loud and clear, with text prompts on a large, bright back-lit display. This makes it ideal if you anticipate use in noisy or poorly lit settings. And at just 2.1 Kg, it's incredibly small and lightweight so it won't weigh you down.

Versatility

The FR2+ is designed to be highly flexible and versatile. You can tailor your FR2+ to match your specific response protocol.

Infant and child defibrillation

The FR2+ can treat patients of any age. When equipped with special infant/child defibrillator pads that reduce the energy of the FR2+'s shock to a more appropriate level, you can safely treat a child or infant in cardiac arrest.

Conscious monitoring**

The reusable 3-leadwire FR2+ ECG assessment module enables the professional responder to proactively use the FR2+ to assess the rhythm of patients who are not in cardiac arrest, but are conscious, breathing, and suspected of being at risk for cardiac distress. This small, lightweight monitoring solution is ideal for in- and out-of-hospital patient transport, bike medics, situations requiring long or difficult travel, fire rigs that may be staffed on occasion by paramedics, clinics and physicians' offices.



Features marked with ** require the M3860A FR2+ model (With ECG and manual override capability)

Ready to use

The HeartStart FR2+ is reliable and virtually maintenance-free. It automatically performs daily, weekly, and monthly self-tests of the battery, electrical components, and subsystems. A highly visible status indicator shows at a glance whether the FR2+ has passed its last self-test and is ready for use.

Reliable analysis technologies

The FR2+ uses Philips' patented SMART Analysis algorithm to assess the patient's heart rhythm and determine if it is shockable. Research shows that for some patients, especially those in long-duration cardiac arrest, CPR prior to defibrillation may provide more benefit. At the discretion of the medical director, the defibrillator's SMART CPR algorithm can be configured to look at characteristics of a shockable heart rhythm and compute the likelihood that spontaneous circulation will return with an initial shock. Depending on the results, the caregiver is advised to provide either an immediate defibrillation shock or CPR first followed by a shock.

Quick Shock

Research also suggests that the opportunity for survival can be improved if interruptions to CPR are minimized. The Quick Shock capability of the FR2+ readies the defibrillator within seconds to deliver a shock, reducing the time between CPR and defibrillation. No other manufacturer can deliver a shock as quickly after CPR.

Proven therapy

When a shock is advised, the FR2+ prompts the user to deliver Philips' low-energy SMART Biphasic waveform, a highly effective defibrillation therapy that is also gentle to the heart. No other external defibrillation therapy has been supported by more published studies.

Comprehensive data management

To maintain high standards of response timeliness and quality, many organizations place a high premium on detailed retrospective event reporting. The FR2+ facilitates this with an optional data card that captures up to 8 hours of patient heart rhythm, defibrillator use and, if configured, voice recording of the event. The removable data card lets you produce a variety of reports using the HeartStart Event Review suite of software products, without taking your FR2+ out of service.

Realistic training


To help responders maintain their skills, Philips gives you a wide range of training solutions. When equipped with training pads and the rechargeable training and administration pack, you can use your FR2+ as a trainer. Its shock delivery capability is disabled while you train with ten scenarios. When you connect the FR2+ to a simulator or a special manikin that simulates an ECG rhythm, you can elicit realistic FR2+ responses to various heart rhythms for even more refined training. Philips also offers an economical AED trainer that comes with ten scenarios. An optional remote control enables instructors to pace the scenario to match their instruction, and to challenge the trainee with surprise scenario changes. Also available is the comprehensive AED Little Anne training system, featuring a manikin that gives the student feedback on pad placement.



Features marked with ** require the M3860A FR2+ model (With ECG and manual override capability)



M3860A WITH ECG and manual override (Display)



M3861A WITHOUT ECG or manual override (Display)

FR2+ Defibrillator

Product Specifications

Defibrillator	
Defibrillator Model	With ECG M3860A With Text Only Display M3861A
How Supplied	Defibrillator; instructions for use, standard battery, defibrillator pads (2 pair), quick reference.
Waveform	Truncated exponential biphasic. Waveform parameters adjusted as a function of patient impedance.
Energy	Single energy output. Adults: 150 J nominal into a 50 ohm load; Infant/Child: 50 J nominal into a 50 ohm load.
Charge Time From Shock Advised	Typically less than 10 seconds.
Charge Time in Manual Mode	Typically less than 5 seconds.
Shock-to-Shock Cycle Time	Typically less than 15 seconds (including analysis time) in automated mode.
Protocol	Text and voice prompts guide user through protocol. Follows pre-configured settings. Can be modified with the M3864A training & administration pack.
Shock Delivery	Via defibrillator pads placed in anterior-anterior (lead II) position for adult defibrillation and anterior-posterior for infant/child defibrillation.
Controls	On/off, shock, screen contrast/option buttons.
Indicators	LCD screen, beeper; audio speakers, status indicator; shock button, connector socket LED.
Advanced Mode	Configurable protocol.

Patient Analysis System	
Patient Analysis	Per protocol, evaluates patient ECG and signal quality to determine if rhythm is shockable, and evaluates connection impedance for proper defibrillation pad contact.
Sensitivity/Specificity	Meets AAMI DF-80 guidelines and AHA recommendations for adult defibrillation (Circulation 1997;95:1677-1682).
SMART CPR (configurable)	Enables support of an automated or user-initiated CPR interval prior to defibrillation. The AUTO 1 and AUTO 2 settings automate the decision of whether to provide CPR first or deliver a shock first based on the amplitude and frequency of the presenting shockable rhythm. Once the decision is made, FR2+ provides the responder with the appropriate prompts. SMART AUTO 1: advises CPR for patients with a presenting rhythm typical of very long-duration cardiac arrest. SMART AUTO 2: advises CPR for an expanded group of patients inclusive of those in Auto 1, having a rhythm typical of long duration cardiac arrest. USER: user initiated CPR Pause interval. Supports a protocol under which the responder decides whether to perform CPR first. A Pause-for-CPR button is enabled, and can be pressed at the responder's discretion.
Quick Shock	Able to deliver a shock in typically less than 10 seconds after the end of a CPR interval.



M3861A WITHOUT ECG or manual override



M3860A WITH ECG and manual override

Physical	
Size	
Height	6.6 cm. (2.6 inches)
Width	21.8 cm. (8.6 inches)
Depth	21.8 cm. (8.6 inches)
Weight	
With battery	2.1 kg. (4.7 pounds)
Without battery	1.8 kg. (3.9 pounds)

ECG Display (M3860A)	
Screen	High-resolution LCD with bright back-light.
Screen Dimensions	2.8 inches wide x 2.3 inches high (7.0 cm x 5.8 cm).
Display Range	Differential: ± 2 mV full scale (nominal).
Sweep Speed	23 mm/second (nominal).
Frequency Response	1 Hz to 20 Hz (-3dB) (nominal).
Sensitivity	1.16 cm/mV (nominal).
Heart Rate	30 to 300 beats per minute updated each analysis period during monitoring.
Monitored Lead	Anterior-anterior (lead II) placement with adult defibrillation pads or ECG Assessment Module (M3860A only).

FR2 Series Standard Battery (M3863A)

Type	12VDC 4.2 Ah lithium manganese. Disposable, recyclable, long-life, primary cells.
Capacity	Minimum 300 shocks or 12 hours operating time.
Install-by Date	Battery is labeled with an install-by date at least 5 years from date of manufacture.
Standby Life	Defines how long the battery will power the AED in standby operation within the standby temperature range (one battery insert test and no uses), 4 years minimum when battery is installed by the install-by date (5 years typical).

Environmental/Physical Requirements

Sealing	Meets IEC529 class IP54 with battery and data card tray installed.
Temperature	Operating: 0° - 50° C (32° - 122° F) Standby: 0° - 43° C (32° - 109° F) Standby applies to AED with battery installed and stored with defibrillation pads.
Humidity	Operating: 0% to 95% relative humidity (non-condensing) Standby: 0% to 75% relative humidity (non-condensing)
Altitude	-500 to 15,000 feet per MIL-STD-810E 500.3 Procedure II.
Aircraft	Method: RTCA/DO-160D: 1997 Section 21 (Category M - Charging).
Shock/Drop Abuse Tolerance	1 meter any edge, corner or surface. MIL-STD-810E 516.4 Procedure IV.
Vibration	MIL-STD-810E 514.4-17
EMI	Requirements: Tested to EN60601-1-2 Limits Radiated: Method EN55011: 1998 Group 1 Level B Immunity: Method EN61000-4-3:1998 Level 2

HeartStart Adult Defibrillation Pads

Configuration	DP2: two pair or DP6: six pair.
How Supplied	Disposable self-adhesive pads with cable and connector.
Surface Area	Meets AAMI DF-80 guidelines.
Cable Length	Approximately 122 cm (48 inches).

HeartStart FR2 Infant/Child Defibrillation Pads (M3870A)

Patient	Under 8 years or 25 kg (55 pounds).
Defibrillator Compatibility	FR2-series (FR2 and FR2+) automated external defibrillator only.
Configuration	M3870A FR2 Infant/Child Reduced Energy Defibrillator Pads. 1 set per package.
How Supplied	Disposable self-adhesive pads with cable and connector.
Energy Delivered	Reduces defibrillator shock energy to nominal 50 Joules into a 50 ohm load.
Surface Area	Meets AAMI DF-80 guidelines.
Cable Length	Approximately 122 cm (48 inches).

Medical Control/Recording Features

Standard Event Review	Elapsed time and number of shocks are displayed on screen.
Enhanced Event Review	Optional data card (M3854A) expands the above on-screen event review capabilities. Review chronological events in detail including ECG. 8 Hours of event & ECG recording or one hour if voice recording is activated.

Automated and User Activated Self-Tests

Automatic Self-Tests	Tests internal circuitry, waveform delivery system, and battery capacity. Verifies calibration of key circuits monthly.
Automatic Self-Test Frequency	Daily when stored within operating environmental conditions.
Status Indication	Dynamic visual and audible indication of self-test results. Indicates device readiness.
Battery Insertion Test	Upon battery insertion, extensive automatic self-tests and user interactive tests check device readiness. Verifies calibration of key circuits.
Automatic Standby Temperature Monitoring	Instrument automatically monitors temperature and warns user if device is stored outside of standby temperature range.

FR2+ Training and Administration Pack (M3864A)

Function	Places FR2+ in scenario-based training mode and disables energy delivery. 10 real-world scripts provided. Permits modification of preprogrammed FR2+ protocol.
Type	12VDC 1.1 Ah Rechargeable Nickel Metal Hydride.
Capacity	Minimum 4 hours training time.
Recharge Time	90 minutes to full capacity using M3855A charger (sold separately).

Its host of features makes the FR2+ an extremely versatile and rugged defibrillator; enabling responders to potentially save the life of a coworker, friend, patient or fellow citizen under demanding circumstances.

Lightweight
Weighs just under 5 pounds fully equipped

Capable
Can be used on anyone of any age, including infants and children

Rugged
Designed for use in extreme environments and conditions

Ever-ready
Powered by either a long-life disposable or a rechargeable battery

Intuitive
1-2-3 operation, with text and voice prompts

Effective
Clinically proven low-energy waveform

HomeHeart Care

Phone: 0845 299 0304
Fax: 0161 601 0204
Email: info@homeheart.co.uk
Web: www.homeheart.co.uk

Home Heart Care Ltd
Sunrise House
Hulley Road
MACCLESFIELD
SK10 2LP
United Kingdom