



PHILIPS

Device / Pads Compatibility

	DP2/DP6	M3713A	M3716A	M3717A	M3718A	M3719A	M3501A	M3502A	M3503A	M3504A	M3870A	M5071A	M5072A	989803139261	Laerdal brand Pads
Pads Device	AED Pads	Adult Pads with Plug Connector	Adult Radiolucent Pads with Plug Connector	Pediatric Pads with Plug Connector	Adult Radio-transparent Pads with Plug Connector	Pediatric Radio-transparent Pads with Plug Connector	Adult Pads with Barrel Connector (AAMI)	Adult Pads with Barrel Connector (IEC)	Pediatric Pads with Barrel Connector (IEC)	Pediatric Pads with Barrel Connector (AAMI)	FR2 Infant/Child Pads	Adult SMART Pads Cartridge	Infant/Child Smart Pads Cartridge	SMART Pads II	Any Laerdal Part Number 902XXX
ForeRunner/FR	X	X	X		X									X*	
FR2/FR2+	X	X	X		X						X			X*	
First Aid/Home/HS1												X	X		
FRx	X*	X*	X*		X*									X	
MRx/XL/XLT/ 4000 with M3507A Barrel Connector							X	X	X	X					
MRx/XL/XLT /4000 with M3508A Plug Connector	X	X	X	X	X	X								X*	
CodeMaster (without 05-10200 adapter)							X	X	X	X					
Any device using 05-10000 (MPC), 05-10100 (ZOLL), or 05-10200 (HP) adapter	X	X	X	X	X	X								X*	
Laerdal HS 3000 with 945005 barrel							X	X	X	X					
Laerdal 3000 with 945003 plug	X	X	X	X	X	X								X*	
Any device using Laerdal manufactured cable 920XXX with Snap Connectors															X

* Pads will operate with device for therapy but cannot be left attached while device is in standby mode or while turned off